
C-REX

EN	Instructions for Use
SV	Bruksanvisning
DA	Brugsanvisning
DE	Bedienungsanleitung

EN

CONTENTS

COMPONENTS.....2
PRODUCT INFORMATION.....5
STEP BY STEP - OVERVIEW.....9
STEP BY STEP.....26
CUSTOMER SERVICE.....51

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




READ THE FOLLOWING INFORMATION CAREFULLY PRIOR TO USE.

This Instructions for Use is designed to provide product information. It is not intended to supersede local institutional protocols or professional clinical judgement concerning patient care.

EN
SV
DA
DE

PRODUCT INFORMATION

INTENDED PURPOSE

C-REX LapAid Cath / C-REX LapAid

C-REX LapAid Cath/C-REX LapAid are sterile, single use, surgical invasive short-term devices that facilitate inverted serosa-to-serosa adaptive anastomosis after resection/resections of colon.

Additionally, the C-REX LapAid Cath variant offers an optional catheter-based anastomotic integrity test, that can provide the surgeon with feedback on intestinal apposition to support further evaluation of the anastomosis.

The devices are suitable for abdominal open surgery, laparoscopic surgery and/or robotic surgery.

The devices can be used for segmental colectomy of descending colon for end-to-end, side-to-side and side-to-end anastomoses.

C-REX RectoAid Cath / C-REX RectoAid

C-REX RectoAid Cath/C-REX RectoAid are sterile, single use, surgical invasive short-term devices that facilitate inverted serosa-to-serosa adaptive anastomosis after resection/resections of colon and/or rectum.

Additionally, the C-REX RectoAid Cath variant offers an optional catheter-based anastomotic integrity test, that can provide the surgeon with feedback on intestinal apposition to support further evaluation of the anastomosis.

The devices are suitable for abdominal open surgery, laparoscopic surgery and/or robotic surgery.

The devices are designed for transanal use for rectal or sigmoid resections to perform end-to-end and side-to-end anastomoses.

C-REX Test Device

The optional accessory C-REX Test device is a sterile, single use, surgical invasive short-term device that during surgery facilitates sizing of the lumen of the intestine, to guide on the proper size of C-REX anastomosis device to use.





DESCRIPTION

Mode of Action: Inverted serosa-to-serosa adaptive anastomosis.

Operating Principle: After resection of colon and/or rectum due to, for instance, tumour, a transient surgical invasive instrument is used to place, in intestine, an anastomotic ring, which is a short-term surgical invasive device, joining colonic ends together. Once Mode of Action is achieved, i.e. the anastomosis is healed at approximately 10 days after surgery, the anastomotic ring loosens from the inside of the intestinal wall and is expelled per rectum.

C-REX LapAid Cath / LapAid and C-REX RectoAid Cath / RectoAid are available in three different sizes and two different variants, i.e. with and without catheters. Please see table on next page.

C-REX Test Device could help for sizing of the intestine. The dimension of C-REX Test Device that fits the intestine gives user the proper size of C-REX device to use, e.g. the first diameter of C-REX Test Device engraved with 26 corresponds to C-REX LapAid Cath 26 / LapAid 26 and C-REX RectoAid Cath 26 / RectoAid 26.

Description	Size	Article No.	Example picture
C-REX Test Device	N/A	A25270	
C-REX LapAid Cath	Size 26 Size 29 Size 32	A26600 A29600 A32600	
C-REX LapAid	Size 26 Size 29 Size 32	A26610 A29610 A32610	
C-REX RectoAid Cath	Size 26 Size 29 Size 32	A26100 A29100 A32100	
C-REX RectoAid	Size 26 Size 29 Size 32	A26110 A29110 A32110	

LapAid – The instrument that is used to apply LapAn in the intestine.

LapAn – A ring-shaped surgically invasive device intended for short-term use, which can be used for both end- and side-anastomosis. An elastic O-ring squeezes the intestinal wall to LapAn.

RectoAid – The instrument that is used to apply RectoAn in the distal intestine.

RectoAn – A ring-shaped surgically invasive device intended for short-term use, which can be applied in distal end-anastomoses. An elastic O-ring squeezes the intestinal wall to RectoAn. RectoAn is placed in the distal part of rectum after resection. The coupling function of RectoAn then makes it possible to perform anastomosis with the other ring, i.e. LapAn, on the proximal intestinal end.

Cath – Through the catheters it is possible to evaluate the anastomosis by conducting anastomotic integrity test.

Integrity test – The integrity pressure test using C-REX catheters is optional and is not intended to supersede clinical judgement concerning any actions to be taken for individual patient care.

Construction of a C-REX anastomosis is such that a negative, traditional air-leak test may give a false negative as the lumen is separated from healing zone. However, it would detect an occult defect outside the anastomosis were there to be an injury to proximal or distal bowel. A positive traditional air-leak test would indicate a failed anastomosis. There is no necessity to undertake an integrity pressure test as the device operates equally with or without catheters. However, integrity pressure testing may provide additional reassurance when constructing the anastomosis confirming good apposition of tissue between proximal and distal bowel.

LIMITATION

C-REX LapAid Cath can only be used when the remaining distal part of the intestine is shorter than the length of the catheters, i.e. 40cm.

C-REX LapAid Cath / LapAid can be used for segmental colectomy on descending colon.

The proximal limit of C-REX RectoAid Cath/RectoAid is <25cm above the anal verge and the distal limit is the anal sphincter.

CONTRAINDICATIONS

- Do not use C-REX device for gastric and oesophageal anastomoses.
- Do not use C-REX device on patients with stenosis or other obstructions distal to the anastomosis.
- Do not use C-REX device for patients who have not undergone bowel preparation using locally agreed purgatives/laxatives prior to surgery.
- Do not use C-REX device where the tissue thickness cannot be comfortably compressed. If the device is used on tissue that cannot be comfortably compressed an inadequate anastomosis could be formed resulting in leakage, inadequate hemostasis or improper healing.

SIDE EFFECTS

- No side effect has been identified during clinical investigations.
- Failure to properly follow **Instructions for Use** may lead to serious surgical consequences, such as bleeding, serosa rift and infection.

CAUTION

C-REX device should only be used according to Instructions for Use.	
C-REX device should only be applied by surgeons whom are adequately trained in using C-REX devices. Inadequate or no training might lead to potential harm.	
Patients shall undergo bowel preparation using locally agreed purgatives/laxatives prior to surgery.	
C-REX device is sterile if the package is dry, unopened and undamaged. The Tyvek lid on the blister is the sterile barrier. Do not use if the package seal is broken or damaged. Do not open the sterile barrier in a non-controlled environment. Microbiologically contaminated medical devices can cause infections in patients.	
C-REX device is intended for single patient use only. Do not resterilize. Dispose opened and unused portions. Attempts to reprocess, resterilize and/or reuse the device can lead to malfunction, transmission of disease and/or negative biocompatible reactions.	 
Dispose C-REX device if mishandling has caused possible damage or contamination or if the device is past its expiration date.	
User has to examine that the anal canal is free from stenosis before choosing this procedure, to ensure that the anastomotic ring can pass through the anal canal.	
During the healing process, the intestinal lumen is limited by the anastomotic ring. To minimise the pressure on the ring, a low residue diet and osmotic laxative is recommended in the early-postoperative period. This can be discontinued when the anastomotic ring has been expelled.	

INFORMATION TO PATIENT

Medical staff has to inform patients on the following:

- To minimise the pressure on the ring, a low residue diet and osmotic laxative is recommended in the early-postoperative period. This can be discontinued when the anastomotic ring has been expelled.
- The external catheters should be kept clean by use of soap and water (when catheter variant is used).
- The patient shall contact medical staff if the anastomotic ring has not left the body within 20 days after the surgery.
- The anastomotic ring shall be disposed as household waste. Avoid flushing down the anastomotic ring in toilet.

DISPOSAL OF SINGLE USE DEVICE

After surgery the device should be disposed according to institutional guidelines for disposal of medical waste.

TECHNICAL PARAMETERS

Device	Size	Overall length (mm)	Outer diameter of Anastomotic ring (mm)	Inner diameter of ring (mm)	Length of Catheters (mm)
DMC Guide handle	-	485	-	-	-
C-REX LapAid Cath / C-REX LapAid	26	363	35	17	440 / -
	29	363	38	20	440 / -
	32	363	41	23	440 / -
C-REX RectoAid Cath / C-REX RectoAid	26	449	35	17	300 / -
	29	449	38	20	300 / -
	32	449	41	23	300 / -

STORAGE

- C-REX device should not be used after expiration date. The duration of validity is indicated with the symbol .
- C-REX device should be stored in a dry and clean indoor environment at normal room temperature.

STERILIZATION

C-REX device has been sterilized using ethylene oxide.

SPECIFICATION OF MATERIALS

ANASTOMOTIC RING
Thermoplastic Polyurethane Elastomer (TPU)
Polyether-etherketone (PEEK)
Silicone rubber
Stainless steel

- The C-REX products are Latex free and free from Phthalates.

LEGEND OF SYMBOLS USED IN INSTRUCTIONS FOR USE AND LABELING

	MEDICAL DEVICE		STERILIZED USING ETHYLENE OXIDE
	CATALOGUE NUMBER		SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING OUTSIDE
	BATCH CODE		SINGLE STERILE BARRIER SYSTEM
	UNIQUE DEVICE IDENTIFICATION		DO NOT USE IF PACKAGE IS DAMAGED
	MANUFACTURER		KEEP DRY For storage.
	DATE OF MANUFACTURE		SINGLE USE ONLY. DO NOT RE-USE
	USE BY DATE. YYYY-MM Do not use any device with an expired date!		DO NOT RESTERILIZE
	CAUTION Consult Instructions for Use for important cautionary information.		The device is CE marked according to Medical Device Directive 93/42/EEC
		0413	


INSTRUCTIONS FOR USE

This Instructions for Use is designed to serve only as a general guideline; it is not a reference to surgical techniques nor intended to supersede local institutional protocols or professional clinical judgement concerning patient care.

This Instructions for Use shows size 32 as examples; the instruction is valid for all sizes.

STEP BY STEP - OVERVIEW

Sizing the intestine using C-REX Test Device, placing the LapAn in the proximal intestine and performing integrity test, are identical for both C-REX LapAid Cath / LapAid and C-REX RectoAid Cath / RectoAid procedures.

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


EN

1 Check the size of intestine (optional)

- a Open Tyvek pouch and Tyvek lid on blister of C-REX device according to hospital routine. C-REX Test Device could help for sizing of the intestine.

 Do not apply excessive force to fit C-REX Test Device into the intestine. Excessive force may lead to ruptured serosa. Always check the integrity of the serosa after use of C-REX Test Device.

- b The dimension of C-REX Test Device that fits the intestine gives user the proper size of C-REX device to use (Size 32 as an example).

 When the diameters of the proximal and distal ends of the intestine differ, use C-REX device with the size that is appropriate for the smaller intestinal diameter.

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

EN

4 Fix anvil in the intestine

The suture shall be located as close as possible to anvil shaft to avoid leakage.

- a Tighten the purse string suture around anvil shaft, the intestine covers round the shaft and the intestinal lumen is closed, so no intestinal content leaks out.

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

EN

5 Connect LapAid to anvil

When anvil with preassembled LapAn is properly fixed with the purse string suture, connect LapAid to anvil shaft. The locking of anvil shaft to LapAid is carried out along anvil connecting pin. Insert anvil connecting pin of LapAid into anvil shaft. User should notice a "click" when anvil snaps over anvil connecting pin.

a End-to-end anastomosis

b Side-to-end anastomosis

User should hold anvil shaft but not the intestine when connecting LapAid to anvil.

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
6 Secure LapAid to the intestine

- a** Rotate adjusting knob clockwise, so LapAid approaches to anvil, which is placed in the intestine (Figure 6a - c).

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

EN

c Side-to-end anastomosis

d The movement can be monitored in the indicator window placed on the upper side of LapAid. The indicator moves into the green range while rotating. Stop rotating and wait about 2 minutes, allowing fluids in tissue to displace. Then continue rotating adjusting knob until it stops.

Check that the indicator on the upper side of LapAid handle is within the green range before firing.

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
EN

7 Fire LapAid

Squeeze trigger handle with firm and steady pressure as far as it will go. This manner will assure that LapAn is fixed to the intestinal wall and detached from anvil, while a circular knife resects excess tissue and cuts the intestine open, making it possible to bring anvil out of the intestine.

User should notice tactile and audible feedback for cutting, which indicates the finish of this step. Release the trigger handle and it will then return back to the original position.

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8 Open up the intestine and remove the anvil

- a Rotate the adjusting knob counter clockwise until it stops, while rotating, the trigger handle goes automatically to the squeezed position (Figure 8a). The attachment of anvil loosens from LapAn and could then be removed. Anvil remains attached to LapAid. Thus LapAn is installed to the proximal intestine, i.e. one half of the anastomosis is created (Figure 8b and 8c).

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
b End-to-end anastomosis

För 8 a, b and c.

If necessary, for instance due to serosa rift, remove the anastomotic ring after placement. The intestine has to be resected again so that unaffected intestine is used with a new C-REX device for anastomosis.

c Side-to-end anastomosis

For C-REX LapAid Cath / LapAid go to Step 9.

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
EN

9 Place 2nd LapAn in the distal intestine

Repeat step 2 to 8, install the 2nd LapAn to the distal intestine.

When step 2 to 9 are completed, the intestinal ends should look like the picture in Figure 9 (if end-to-end anastomosis is performed).

10 Bring DMC Guide handle/DM

- a Once the two LapAn rings are fixed in the intestine, they should be connected together by using the coupling segment DMC (part of DMC Guide handle, for C-REX LapAid Cath) or DM (for C-REX LapAid).

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
EN

b Coupling segment DM

11 Please follow either step 11a-c (for C-REX LapAid Cath) or 11d (for C-REX LapAid)

- a** The guide shall be led into the distal intestine and out through anus. It is recommended to insert one finger through the anus to lead top cone out from anus.

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
EN

b Snap DMC and LapAn together.

c DMC Guide handle is pulled out through anus.

Check the attachment of DMC by pressing DMC to LapAn carefully and evenly. Make sure that the rim of DMC is parallel to the rim of LapAn.

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
EN

d Fix DM to distal LapAn

Connect DM by inserting DM into LapAn, and snapping DM and LapAn together. Check attachment of DM by pressing DM to LapAn carefully and evenly. Make sure that the rim of DM is parallel to the rim of LapAn.

Avoid unnecessary stretch of the distal intestinal end to prevent tearing.

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12 Connect the proximal and distal LapAn

Carefully inspect the proximal and distal intestine before connection to ensure that the mesentery is aligned without twisting.

- a** The proximal and distal intestine installed with LapAn rings and DMC/DM are shown.

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

b Side-to-end anastomosis

- c The proximal and distal LapAn rings are connected with a simple locking procedure as follows:

1. Place thumbs 180° from each other on one side of anastomosis.
2. Fingers should be parallel to thumbs and placed on the opposite LapAn.
3. Apply even pressure on both LapAn rings until they are locked together.
4. Move fingers 90° and squeeze again to ensure the locking.

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
d End-to end anastomosis

For 12 a,b,c and d.

Upward traction on the proximal intestinal end should be avoided to prevent tearing. Loosen the left flexure if needed.

For 12 a,b,c and d.

It is important that there is a serosa-to-serosa contact circumferentially. Check that LapAn rings are properly locked to each other. Figure 12d show a section view of the completed anastomosis (with catheters).

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For 12 a,b,c and d.

Do not put any stitches around the anastomosis.

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Für 8 a, b, c und d.

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EN

**Go to Step 19 for C-REX LapAid Cath;
Anastomotic integrity pressure test
or to Step 20 for C-REX LapAid.
See page 50.**

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

EN

13 Connect RectoAn Anvil or RectoAn Anvil Cath to proximal anastomotic ring

- a Take out RectoAn Anvil or RectoAn Anvil Cath from the package.
Insert RectoAn Anvil to the proximal LapAn ring.

- b LapAn and RectoAn are connected with a simple locking procedure as follows:

1. Place thumbs 180° from each other on the proximal anastomotic ring.
2. Fingers should be parallel to thumbs and placed on RectoAn.
3. Apply even pressure on both rings until they are locked together.
4. Move fingers 90° and squeeze again to ensure the locking.

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
EN

- c When RectoAn Anvil is placed correctly in the proximal end, intestinal content cannot flow into abdominal cavity (waterproof). Place then the intestine with the anvil close to the distal intestine end.
At robotic/laparoscopic surgery, the pfannenstiel-incision shall be closed and insufflation shall be started.

14 Insert RectoAid through anus

Place the patient in Lithotomy position. It is recommended to do a rectal examination at this time. Lubricate the anus with lubricant and insert the outer tube of RectoAid through anal canal.

Insert RectoAid carefully to avoid damage of anal canal and rectum.

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
EN

15 Connect RectoAid to RectoAn Anvil

- a When the outer tube of RectoAid has passed the anal canal, carefully insert it further up in the rectal stump (distal intestine). Push RectoAid forward until the distal intestine is stretched out over RectoAid outer tube.

Violent pushing of RectoAid shall be avoided when inserting RectoAid up.

Upward traction on the intestinal ends should be avoided to prevent tearing.

- b Rotate adjusting knob counter clockwise until it stops, so anvil connecting pin of RectoAid protrudes and penetrates the closed distal intestine (Figure 15b and 15c).

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EN

c

- d The locking of RectoAid to RectoAn anvil shaft is carried out along anvil connecting pin. Insert anvil connecting pin into anvil shaft. User should notice a "click" when anvil snaps over anvil connecting pin.

User should hold anvil shaft but not the intestine when connecting RectoAid to anvil.

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

EN

16 Secure RectoAn to the intestine

- a** Rotate adjusting knob clockwise, so RectoAid approaches to RectoAn Anvil, which is placed on the LapAn on the proximal intestine (Figure 16a and 16b).

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

b

- c The movement can be monitored in the indicator window placed on the upper side of RectoAid. The indicator moves into the green range while rotating.

- d Stop rotating and wait about 2 minutes, allowing fluids in tissue to displace. Then continue rotating adjusting knob until it stops. Now the proximal and distal intestine adapt to each other.

Check that the indicator on the upper side of C-REX RectoAid is within the green range before firing. Carefully inspect the proximal and distal intestine before firing to ensure that the mesentery is aligned without twisting.

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EN

17 Fire RectoAid

Squeeze trigger handle with firm and steady pressure as far as it will go. RectoAn is fixed to the intestinal wall and automatically detached from anvil, while a circular knife resects excess tissue and cuts the intestine open, making it possible to bring RectoAn anvil out. User should notice tactile and audible feedback for cutting, which indicates the finish of this step. Release the trigger handle and it will then return back to the original position.

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
18 Remove the RectoAn Anvil and pull out RectoAid

- a The anvil remains attached to RectoAid. Pull out RectoAid through anus. Thus RectoAn is installed to the distal intestine, and in the meanwhile it connects to the proximal anastomotic ring, i.e. the anastomosis is created. Figure 18a shows the section view of the completed anastomosis (with catheters).

It is important that there is a serosa-to-serosa contact circumferentially.

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Step 19. Anastomotic integrity pressure test

(for C-REX LapAid Cath and C-REX RectoAid Cath only)

To confirm that the two anastomotic rings are properly locked to each other, check the anastomotic integrity pressure between the intestinal ends by using the catheters.

PRESSURE GAUGE: Use an ordinary digital pressure gauge with a range of 0 - 300 mbar and accuracy ± 10 mbar.

AIR PRESSURE SOURCE: Flow 0.5 ml/sec, max air pressure 300 mbar (example: syringe).

FILTER: Sterile filter (luer lock).

TUBING: Syringe tubing with luer lock end (included in this C-REX device).

Connect the catheter tubing to the syringe tubing. Connect the air pressure source, the pressure gauge and a filter to the luer lock end of syringe tubing via a multi-way stopcock according to Figure 19 above.

Pour sterile saline into the open abdominal cavity so that the saline covers the anastomosis.

Start the measuring function on the pressure gauge, and carefully increase the air pressure (approx. 0.5 ml/sec) through the air pressure source. When air bubbles appear from the anastomosis the integrity pressure is reached. Read the pressure shown on the pressure gauge.

Based on observational data, an integrity pressure above approximately 80 mbar has been found to indicate good apposition of the intestinal ends, however, it remains at the surgeon's discretion to determine whether a lower value is acceptable in the context of an otherwise healthy-looking anastomosis. The integrity pressure test is not intended to supersede clinical judgement concerning any actions to be taken for individual patient care.

Disconnect the pressure gauge, syringe, filter and syringe tubing.

Use surgical aspirator to drain the abdominal cavity from saline.

Step 20. Close the abdomen

Close the abdomen according to local institutional surgical procedure.

Step 21. Measure the length of catheters outside anus (for C-REX LapAid Cath and C-REX RectoAid Cath only)

For patients' convenience, cut the 4-to-1 connector and catheters after surgery, and leave approximately 10 cm long ends from the anal rim.

Measure the distance from the anus to the free ends of the catheters (outside anus) and note it in the patient record.

Monitoring this distance provides a simple way to monitor the anastomotic ring. An increase of the length of the catheters indicates that the anastomotic ring has loosened. After approximately 10 days the anastomotic ring loosens from the inside of intestine and follow the faeces out.

Patients are not allowed to pull the catheters since it could cause intestinal leakage. On mentally unstable patients, the catheters should be cut as short as they do not protrude outside anus. The cutting procedure may only be done by trained surgeons.

Patients have to contact medical staff if the anastomotic ring has not left the body within 20 days after surgery. Intervention could be necessary to remove anastomotic ring.

Step 22. Postoperative radiological examination of the anastomosis (if necessary)

Please note that the radiological examination of the anastomosis can be performed during operation as well with same protocol as described below.

X-ray examination via the catheters could be done within the first 48 h after surgery but only when anastomotic leakage is suspected.

PRESSURE GAUGE: Use an ordinary digital pressure gauge with a range of 0-300 mbar and accuracy ± 10 mbar.

AIR PRESSURE SOURCE: Max air pressure 100 mbar (example: syringe, ball pump).

NEEDLE: Sterile needle 21G.

FILTER: Sterile filter.

CLAMP: Artery forceps.

TUBING: Sterile plastic tube with luer lock ends. Inner diameter 1 - 6 mm.

X-RAY EQUIPMENT: for instance a C-arm X-ray.

WATER SOLUBLE CONTRAST MEDIA: for instance an iodinated, soluble X-ray contrast.

Before X-ray examination:

Cut off the 4-to-1 connector on the four catheters if it is not done before.

Make sure that the catheters are clean:

Dip three of the catheters into sterile saline and aspirate (suck) with a syringe in the fourth catheter.

Repeat the procedure in all catheters.

During X-ray examination:

Dip three of the catheters into water-soluble contrast media and aspirate (suck) with a syringe in the fourth catheter. Contrast media is thus aspirated into the catheters to the spatium between the intestinal ends.

Clamp the three catheters dipped in contrast media.

Connect the fourth catheter to the pressure gauge.

Raise the pressure in the spatium between the intestinal ends by infusing air in the fourth catheter, while reading the pressure showed on the pressure gauge.

In the meanwhile conduct X-ray examination.

The maximum pressure should not exceed half of the integrity pressure measured during surgery in step 19, but max 100mbar.

In case of anastomotic leakage, a leakage of contrast media could be observed under X-ray examination, and no pressure could be detected.

After X-ray examination:

Dip three of the catheters into sterile saline and aspirate (suck) with a syringe in the fourth catheter to clean the catheters from contrast media.

Repeat the procedure in all catheters.

This examination only detects leakage in the anastomosis.

The contrast media should be applied only by vacuum.

It is easy to break up the healed area in the anastomosis, if this step is not carried out carefully.

Step 23. Evacuation of the anastomotic ring

The anastomotic ring is expelled per rectum at approximately 10 days after surgery. The anastomosis is considered healed at this stage.

Post-operative patient management (such as drainage, diet, and patient mobilization) may affect ring expulsion time to be shorter or longer than 10 days. Post-operative liquid diet may delay the time for ring evacuation.

Patients have to contact medical staff if the anastomotic ring has not left the body within 20 days after surgery. Intervention could be necessary to remove the anastomotic ring.

REPORTING

Any complaint or serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user is established.

CUSTOMER SERVICE

If you have any questions, please contact CarpoNovum, address and phone number found below.

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CarpoNovum AB
CarpoNovum AB | Olofsdalsvägen 10 | SE-302 41 Halmstad | SWEDEN